

Leaders in
Diagnostic Solutions

Audiology Equipment





• Amplivox 240 Audiometer

On Tender!

RT274-01-022

Included:

Carry bag and batteries

Technical Specs



- Compact and portable
- Air conduction (125 – 8,000 Hz)
- Bone conduction (250 – 8,000 Hz)
- Masking
- Optional battery power function
- Automatic air conduction functionality
- Talk over function
- Results transferable to PC



• Auditdata Measure

On Tender!

RT274-01-064

Technical Specs



- Real Ear Measurement/Verification equipment
- On ear and test box speech mapping
- RECD
- HIT (Hearing Instrument Testing)
- Variety of stimuli, including ISTS
- Workflows customisable for hospital protocol
- PC-based
- Optional add-ons:
 - Two-channel diagnostic, PC-based audiometry module

audiostar^{pro}TM

gsi
Grason-Stadler

By the no.1 trusted brand



Product Features

- ACT test (non-language specific speech in noise test) unique feature!
- True 2-channel clinical audiometer
- The preferred choice for paediatric assessments
- Air, bone, sound field
- Speech (live and pre-recorded)
- Special tests incl. TEN to identify dead regions
- High frequency audiometry (up to 20 kHz) for ototoxicity monitoring
- Results transferable to PC (PC-controlled license available as add-on)

R165'000*

incl. VAT

SAVE R28'000!

*Prices are subject to exchange rate fluctuations and may change without prior notice

Specs & Motivation



pello™

By the no.1 trusted brand

gsi
Grason-Stadler



Product Features

- ACT test (non-language specific speech in noise test)
- Air, bone, masking
- Speech (live and pre-recorded)
- Results transferable to PC
- Special tests incl. TEN to identify dead regions
- High frequency audiometry (up to 20 kHz) for ototoxicity screening

unique feature!

R99'000*

incl. VAT

SAVE R13'000!

*Prices are subject to exchange rate fluctuations and may change without prior notice

Specs & Motivation



Auditdata

Measure

The next generation of portable audiometry



Product Features

- Ambient Noise Microphone
- 2 channel diagnostic audiometer
- PC-based
- Portable
- Air, bone, speech (live and pre-recorded)
- Sound field enabled
- Special tests available
- High frequency audiometry (up to 16 kHz) for ototoxicity monitoring

unique feature!

R147'500*
incl. VAT

SAVE R43'000!

*Prices are subject to exchange rate fluctuations and may change without prior notice

Specs & Motivation



tymptestar^{pro}

gsi
Grason-Stadler

The gold standard for middle ear assessments



Product Features

- GSI Durability
- Wideband Tympanometry unique feature!
- 226, 678 and 1,000 Hz probe tones
- Pulsed acoustic reflexes for increased accuracy
- Auto Sequence
- 3 Eustachian Tube Tests
- Reflex Decay and Reflex Latency
- Touchscreen
- Results transferable to PC

R156'000*
incl. VAT

SAVE R20'000!

*Prices are subject to exchange rate fluctuations and may change without prior notice

Specs & Motivation



*Prices are subject to exchange rate fluctuations and may change without prior notice

novus™

gsi
Grason-Stadler

Setting the standard for early detection



Product Features

- CE-Chirp for faster ABR unique feature!
- Objective and accurate
- Ideal for newborn screening
- Flexible, single-device workflow for screening OAE and ABR
- Touch screen display
- Compact hardware design
- Test ABR in both ears simultaneously
- Wireless charging
- Results transferable to PC

R112'000*
incl. VAT

SAVE R15'000!

*Prices are subject to exchange rate fluctuations and may change without prior notice

Specs & Motivation



corti™

gsi
Grason-Stadler

Setting the standard for early detection



Product Features

- Unique patented noise rejection algorithm – ideal for variable noise levels unique feature!
- Diagnostic DPOAE and TEOAE
- Screening DPOAE and TEOAE included!
- Robust design that can withstand the demands of high-volume clinics
- DPOAE up to 12 kHz for ototoxicity monitoring
- Results transferable to PC

R129'000*
incl. VAT

SAVE R17'000!

*Prices are subject to exchange rate fluctuations and may change without prior notice

Specs & Motivation



New Equipment

Booths

Calibrations

Repairs

Consumables

Training

Support



The Importance of Annual Calibration

Your key to accurate diagnosis



Ensure **accurate results** in line with HPCSA guidelines



Keeps equipment **running smoothly** and reduces unexpected issues



Prevents errors and **ensures compliance** with SANS/IEC standards



Extends equipment lifespan and helps avoid costly repairs



Equipment Quotations & Calibrations:

audioequip@batho.co.za

012 763 6816

079 970 5021 (WhatsApp)

Request a quote:

